



Infectious Disease by Nucleic Acid Testing

Specimen requirement

TEST	SPECIMEN	CONTAINER	STORAGE TEMP	SHIPPING TEMP
RESPIRATORY PANEL	SWAB	VIRAL TRANSPORT MEDIA (VTM)	AMBIENT	AMBIENT
GASTROINTESTINAL PANEL	STOOL	STERILE CUP	REFRIGERATE 2-6 C	AMBIENT
SEXUALY TRANSMITTED DISEASE	URINE SWAB	URINE CUP OR VAGINAL SWAB	AMBIENT	AMBIENT
BLISTER PANEL	SWAB	VIRAL TRANSPORT MEDIA (VTM)	AMBIENT	AMBIENT
BLOOD VIRAL PANEL	BLOOD	LAVENDER TOP (EDTA TUBE)	AMBIENT	AMBIENT
MRSA/MSSA	SWAB	VIRAL TRANSPORT MEDIA (VTM)	AMBIENT	AMBIENT

Fertility & Cryogenics Lab
 8635 Lemont Road
 Downers Grove, IL 60516-4805
 Voice (630) 427-0300 Fax (630) 427-0302